

Y-Knot® PRO

Cleatless All-Suture Anchor System

Available with CONMED's Y-Knot® RC and Y-Knot® Flex all-suture anchors, Y-Knot® PRO anchors are designed to improve efficiency by eliminating the need to uncleat the sutures after anchor insertion.



Features and Benefits

Cleatless Delivery

Providing a faster delivery, Y-Knot PRO anchors eliminate the need to uncleat sutures.

Small, Strong All-Suture Anchors

Y-Knot PRO combines an enhanced delivery system with all the benefits of the Y-Knot RC & Y-Knot Flex.

- Small: Designed to improve placement options in limited anatomy.
- Soft: No hard, loose bodies left behind.
- Simple: Y-Knot Flex anchors feature curved, flexible instrumentation while Y-Knot RC Anchors are self-punching.
- Strong: Y-Knot Flex and Y-Knot RC anchors feature pull-out strength of 477N¹ and 602N¹ respectively.

Your Choice of Suture, Ribbon & Tape

Y-Knot PRO anchors are available single, double or triple-loaded with CONMED's #2 Hi-Fi® Suture, Hi-Fi® Ribbon or Hi-Fi® Tape.

Multiple Procedures

Y-Knot PRO Flex anchors can be used for Glenoid instability, PASTA repair, Distal and Proximal Biceps repair, and Lateral Ankle Instability procedures.

Y-Knot PRO RC anchors can be used for Rotator Cuff repair, Massive Rotator Cuff Revision/Augmentation, Superior Capsular reconstruction, Remplissage, Distal or Proximal Biceps repair and Achilles procedures.

¹Barber A., et al. All-Suture Anchors: Biomechanical Analysis of Pullout Strength, Displacement, and Failure Mode. Arthroscopy 2017; 33:1113-1121.




To simplify your all-suture fixation, visit [CONMED.com/Y-KnotPRO](https://www.conmed.com/Y-KnotPRO) or call **1-800-237-0169** to schedule a demonstration.



525 French Road, Utica, NY 13502

Y-KNOT® PRO
CLEATLESS
ALL-SUTURE
ANCHOR SYSTEM



Y-Knot® PRO Flex All-Suture Anchor
W/ TWO 1.3MM HI-FI® RIBBONS

Y-Knot® PRO RC All-Suture Anchor
W/ TWO #2 HI-FI® SUTURES

Ordering Information

Description	Size	Catalog Number
NEW Y-Knot® PRO Flex All-Suture Anchor		
Y-Knot® PRO Flex w/ one #2 Hi-Fi® Suture (White/Black)	1.3mm	YP1301W
Y-Knot® PRO Flex w/ one #2 Hi-Fi® Suture (Blue)	1.3mm	YP1301B
Y-Knot® PRO Flex w/ two #2 Hi-Fi® Sutures (Blue and White/Black)	1.8mm	YP1802
Y-Knot® PRO Flex w/ one 1.3mm Hi-Fi® Ribbon (White/Black)	1.3mm	YP1301RW
Y-Knot® PRO Flex w/ one 1.3mm Hi-Fi® Ribbon (White/Blue)	1.3mm	YP1301RB
Y-Knot® PRO Flex w/ two 1.3mm Hi-Fi® Ribbons (White/Blue and White/Black)	1.8mm	YP1802R
NEW Y-Knot® PRO RC All-Suture Anchor		
Y-Knot® PRO RC w/ two #2 Hi-Fi® Sutures (Blue, White/Black)	2.8mm	YPRC02
Y-Knot® PRO RC w/ three #2 Hi-Fi® Sutures (Blue, White/Black and White/Blue)	2.8mm	YPRC03
Y-Knot® PRO RC w/ one 2mm Hi-Fi® Tape (White/Black), one #2 Hi-Fi® Suture (Blue)	2.8mm	YPRCTW
Y-Knot® PRO RC w/ one 2mm Hi-Fi® Tape (Blue), one #2 Hi-Fi® Suture (White/Blue/Black)	2.8mm	YPRCTB
Y-Knot® PRO RC w/ two 1.3mm Hi-Fi® Ribbons (White/Blue and White/Black)	2.8mm	YPRC02R
Y-Knot® PRO RC w/ three 1.3mm Hi-Fi® Ribbons (Blue, White/Black and White/Blue)	2.8mm	YPRC03R

Local 727-392-6464
Toll Free 1-800-237-0169

CONMED.com
customer_service@conmed.com