




Gastroenterology and Pulmonary Product Catalog




DURA CLIP™
HEMOSTATIS CLIP



BEAMER® ARGON
PROBE



ELIMINATOR® PET
BALLOON DILATOR



BEAMER® ARGON
CE600



To order these and other
innovative products,
call **866-4CONMED** or
visit **CONMED.com**.

TABLE OF CONTENTS

ORDERING INFORMATION

Ordering Information	1
----------------------	---

BILIARY

Papillotomes

TruPass® SW Biliary Papillotome	2
Apollo 3AC® Advanced Cannulation Triple Lumen Papillotomes	3
Apollo 3® Tapered Tip Triple Lumen Papillotomes	3

Cannulae

ProForma® HF 4.5 Cannula	4
--------------------------	---

Guidewires

ACE™ Biliary Guidewire	5
FXWire® Advanced Measurement Guidewire	5
XWire® Guidewire	6

Dilation Balloons

Eliminator® PET Biliary Balloon Dilators	7
--	---

Stone Removal Balloons

FIT® Biliary Stone Balloon	8
DuraGlide™ Triple Lumen Stone Removal Balloons	8

Metal Stents

GORE® VIABIL® SW Biliary Endoprosthesis	9
GORE® VIABIL® LW Biliary Endoprosthesis	10
Flexxus® Endoscopic Biliary Stent	11

ENDOSCOPIC ULTRASOUND

Fine Needle Aspiration

ClearView™ Endoscopic Ultrasound Needle	12
---	----

DILATION

PET Balloon Dilators

Eliminator® PET Balloon Dilators	13
CONMED Balloon Inflation Device	14
Multi-Stage Inflation Pressures	14

Reusable Polyvinyl Dilators

American Dilation System	15
CleanGuide™ Disposable Spring Tip Guidewire	16

HEMOSTASIS

Hemostasis Clips

DuraClip™ Repositionable Hemostasis Clip	17
--	----

Hemostasis Probes

BiCap® Superconductor™ Multi-electrode Bipolar Probes	17
Beamer® Argon Snare Probe	18
Beamer® Argon Probe	18
Beamer® Argon SideFire Probe	18
ABCFlex™ Argon Probes	19
Coaxial Bipolar Adaptors	19
Active Cord	19

Disposable Injection Needles

Click-Tip™ Injection Needle	20
Flexitip™ Disposable Sclerotherapy Needles	20

Ligation Devices

Stiegmann-Goff™ Endoscopic Ligator Kit	21
CONMED Endoscopic Overtube	21
Stiegmann-Goff™ Bandito™ Endoscopic Ligator	22
SpaceBander® Hemorroid Band Ligator	22

BIOPSY

Disposable Biopsy Forceps

Precisor® EXL™ Coated Disposable Biopsy Forceps	23
Optibite® Disposable Biopsy Forceps	23
Pediatric GI	24
Precisor® Hot Disposable Biopsy Forceps	24
Precisor® Broncho® Disposable Biopsy Forceps	24

Disposable Percutaneous Biopsy Instrument

Monotyp® Percutaneous Biopsy Instrument Kits	25
--	----

Disposable GI Cytology Brushes

CONMED Cytology Brushes	26
Cytology Brushes	26

GASTROSTOMY

Gastrostomy Feeding Systems

EnTake® PEG System (No Drugs)	27
EnTake® PEG System (No Drugs): Accessories	27
EnTake® PEG System	28

POLYPECTOMY

Disposable Snares

Beamer® Argon Snare Probe	29
Singular™ Polypectomy Snares	30
Optimizer® Polypectomy Snares	31
Orbit-Snare® Rotatable Polypectomy Snares	32
CONMED Polypectomy Snares	32
Active Cord	32

Specimen Devices

Nakao Spider-Net® Retrieval Devices	33
Optimizer® Multi-Chamber Polyp Trap	33

INFECTION PREVENTION

Biopsy Valves / Introducer

CleanGuard™ Biopsy Valves	34
Biopsy Channel Accessory Introducer	34

Bite Blocks

Scope Saver® Bite Blocks	34
CONMED Bite Blocks	34

Endoscope Valve-Tracking and Pre-Cleaning System

CleanCage™ Endoscope Valve-Tracking System	35
CleanTray™ Endoscope Pre-Cleaner	35

CleanGuard™ Disposable Endoscope Valves

CleanGuard™ Disposable Endoscope Valves	36
---	----

Cleaning Brushes

Channel Master™ Endoscope Cleaning Brushes	37
Whistle® Endoscope Cleaning Brush	37
Blue Bullet® Disposable Cleaning Brushes	38
Reusable Cleaning Brush	38

PULMONARY

Transbronchial Aspiration Needles

WANG® Cytology Needles	39
WANG® Histology Needles	39

Biopsy Forceps

Precisor Broncho® Disposable Biopsy Forceps	40
---	----

Brushes

Cytology Brushes	41
Needle Tip and Peripheral Cytology Brushes	42
Microbiology Brush	42

Cleaning Brushes

Blue Bullet® Disposable Cleaning Brush	43
Channel Master™ Endoscope Cleaning Brush	43
Reusable Cleaning Brush	43

GENERATORS

Generators

Beamer® System CE600	44
BiCap® III Electrosurgical Generator	45
Irrigation Pump	45

Argon Probes

Beamer® Argon Snare Probe	46
Beamer® Argon Probe	46
Beamer® Argon SideFire Probe	46
ABCFlex® Argon Probes	46

LATEX PRODUCTS 47

CATALOG NUMBER INDEX 48

To Place Your Order

Please provide a purchase order number, billing and shipping instructions, product catalog number, description and pricing to:

Customer Service if calling from within the U.S.

Telephone: (866) 4CONMED
Fax: (800) 438-3051
E-mail: customerservice@CONMED.com

International Customer Service

Telephone: +1(727) 214-3000
Fax: +1(315) 735-6235

Terms and Conditions

All prices and specifications are subject to change without notice. All prices are F.O.B. Utica, New York and Lithia Springs, Georgia, (F.O.B. Drop Shipment Points where applicable). Purchase orders should be addressed to CONMED Corporation, 525 French Road, Utica, NY 13502.

CREDIT TERMS: Net 30 days.

Return Goods Policy

All CONMED products, except products that are obsolete, discontinued, or damaged and specialty products, may be returned within 60 days after date of purchase for credit subject to the following:

1. A copy of the invoice number and date of purchase and a Return Authorization Number must accompany the returned goods with the reason for the return stated in writing. Contact CONMED Customer Service for the Return Authorization Number at 1-866-4CONMED.
2. All returned goods would be subject to a minimum of 15% restocking charge. The exact amount of credit to be given, if any, will be determined by a Quality Assurance inspection. Sterile packaged products may be returned for credit only if returned in the original, unopened shelf package.
3. If the return is due to our error in shipping, full credit will be given.
4. All returns should be addressed as indicated by Return Authorization.

Scheduled Shipments

To help maximize the efficiency of your ordering process and inventory control, we offer scheduled shipments.

The advantages include:

- regular delivery schedule
- one purchase order per year
- favorable pricing
- price protection
- free freight for Federal Express Economy Delivery

Quality Policy

At CONMED our passion and responsibility is for the success of our customers and the patients for whom they care.

We empower our customers to succeed by our commitment to integrity, quality, responsiveness and a relentless focus to deliver accessible CONMED solutions. Our talented employees and culture are the foundation of our success.

Papillotomes



TruPass® Short Wire Triple Lumen Papillotome

Sterile, Disposable

Catalog Number	Type	Tip		Cutting Wire		Catheter		Guidewire	Unit
		Length	O.D.	Length	Type	Length	O.D.	Compatibility	
T2043F	Short Nose	8mm	4.3 F	20mm	Monofilament	190cm	7-4.3 F	.035"/.89mm	1/box
T3043F	Short Nose	8mm	4.3 F	30mm	Monofilament	190cm	7-4.3 F	.035"/.89mm	1/box

- Guaranteed Cutting Wire Orientation between 11:00 - 1:00
- 4.3 F .035" Compatible Triple Lumen Tapered Tip
- Smooth, Atraumatic Radiopaque Tip
- Integrated Lumen Technology for Optimal Guidewire Positioning
- Easily identifiable green guidewire exit port
- Exit Marker Band
- Cutting Wire Markers

Papillotomes




Apollo 3AC® Advanced Cannulation Triple Lumen Papillotomes

Sterile, Disposable

Catalog Number	Type	Tip		Cutting Wire		Catheter		Guidewire Compatibility	Unit
		Length	O.D.	Length	Type	Length	O.D.		
7104AC	Short Nose	5mm	4.5 F	20mm	Monofilament	190cm	7-4.5 F	.035"/.89mm	1/box
7105AC	Short Nose	5mm	4.5 F	30mm	Monofilament	190cm	7-4.5 F	.035"/.89mm	1/box

- Guaranteed Cutting Wire Orientation
- 4.5 F .035" Compatible Triple Lumen Tapered Tip
- Color Coded Guidewire Port
- Exit Marker Band
- Cutting Wire Markers



Apollo 3® Tapered Tip Lumen Papillotomes

Sterile, Disposable

Catalog Number	Type	O.D.	Cutting Wire		Catheter		Guidewire Compatibility	Unit
			Length	Type	Length	O.D.		
007104	Short Nose	5.0 F	20mm	Monofilament	190cm	7-5.0 F	.035"/.89mm	1/box
007105	Short Nose	5.0 F	30mm	Monofilament	190cm	7-5.0 F	.035"/.89mm	1/box

■ Cannulae




ProForma® HF 4.5 Cannula

Sterile, Disposable

Catalog Number	Tip	Tip O.D.	Catheter O.D.	Length	Guidewire Compatibility	Unit
Precurved						
050907	Curved Taper Tip	4.5 F/1.5mm	7-5 F/1.7mm	190cm	.035"/.89mm	5/box
050906	Straight Taper Tip	4.5 F/1.5mm	7-5 F/1.7mm	190cm	.035"/.89mm	5/box


Guidewires



ACE™ .025 Biliary Guidewire

Sterile, Disposable

Catalog Number	Description	Guidewire O.D. Inches	Tip Length	Wire Length	Tip Type	Torque Vise	Unit
Straight Guidewires							
A25261S	ACE™ .025 Biliary Guidewire	.025"	10cm	260cm	Straight Tip	Yes	2/box
A25481S	ACE™ .025 Biliary Guidewire	.025"	10cm	480cm	Straight Tip	Yes	2/box
Angled Guidewires							
A25261A	ACE™ .025 Biliary Guidewire	.025"	10cm	260cm	Torquable Angled Tip	Yes	2/box
A25481A	ACE™ .025 Biliary Guidewire	.025"	10cm	480cm	Torquable Angled Tip	Yes	2/box




FXWire® Advanced Measurement Guidewire

Sterile, Disposable

Catalog Number	Description	Guidewire O.D. Inches	Tip Length	Wire Length	Tip Type	Torque Vise	Unit
FX35450S	Standard Wire	.035"	5cm	450cm	Straight Tip	No	3/box
FX35450A	Standard Wire	.035"	5cm	450cm	Torquable Angled Tip	Yes	3/box
FX35260S	Standard Wire	.035"	5cm	260cm	Straight Tip	No	3/box
FX35260A	Standard Wire	.035"	5cm	260cm	Torquable Angled Tip	Yes	3/box

■ Guidewires




XWire® Guidewire

Sterile, Disposable

Catalog Number	Description	Guidewire O.D. Inches	Tip Length	Wire Length	Tip Type	Torque Vise	Unit
035450S	Standard Wire	.035"	5cm	450cm	Straight Tip	No	3/box
035450A	Standard Wire	.035"	5cm	450cm	Torquable Angled Tip	Yes	3/box
035260S	Standard Wire	.035"	5cm	260cm	Straight Tip	No	3/box
035260A	Standard Wire	.035"	5cm	260cm	Torquable Angled Tip	Yes	3/box
035450SS	Stiff Wire	.035"	5cm	450cm	Straight Tip	No	3/box
035450AS	Stiff Wire	.035"	5cm	450cm	Torquable Angled Tip	Yes	3/box
035260SS	Stiff Wire	.035"	5cm	260cm	Straight Tip	No	3/box
035260AS	Stiff Wire	.035"	5cm	260cm	Torquable Angled Tip	Yes	3/box
025450S	Standard Wire	.025"	5cm	450cm	Straight Tip	No	3/box
025450A	Standard Wire	.025"	5cm	450cm	Torquable Angled Tip	Yes	3/box
025260S	Standard Wire	.025"	5cm	260cm	Straight Tip	No	3/box

■ Dilation Balloons



Eliminator® PET Biliary Balloon Dilators

Sterile, Disposable

Catalog Number	Description	Balloon Inflation		Catheter		Guidewire	Unit
		Pressure	Length	O.D.	Length	Compatibility	
000347	4mm x 2cm	150 psi/10.2 atm	2cm	5 F/1.7mm	180cm	.035"/.89mm	1/box
000348	6mm x 2cm	150 psi/10.2 atm	2cm	5 F/1.7mm	180cm	.035"/.89mm	1/box
000349	8mm x 3cm	120 psi/8.2 atm	3cm	5 F/1.7mm	180cm	.035"/.89mm	1/box

Stone Removal Balloons




- One Balloon for Any Size Duct
 - Expandable from 9mm-22mm
- Designed to Reduce Resistance and Bursting
- Dual-Injection Ports
- Choice of proximal or distal injection
- Unlike current market-leading biliary stone balloons, FIT® is not made with natural rubber latex.
- CONMED Biliary Assurance Program

FIT® Biliary Stone Balloon

Sterile, Disposable

Not made with natural rubber latex.

Catalog Number	Balloon		Catheter		Guidewire	Unit
	Size	Max Volume	O.D.	Length	Compatibility	
Standard						
F922SW	9-22mm	4.5mL	7 F/2.3mm	210cm	.035"/.89mm	1/box
F922LW	9-22mm	4.5mL	7 F/2.3mm	210cm	.035"/.89mm	1/box




DuraGlide™ Triple Lumen Stone Removal Balloons

Latex

Sterile, Disposable

	Balloon		Catheter		Guidewire	
Catalog Number	O.D.	Volume	O.D.	Length	Compatability	Unit
Injection Above Balloon						
007123	8.5mm	.75 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box
007124	11.5mm	1.5 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box
007125	15.0mm	3.5 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box
Injection Below Balloon						
007143	8.5mm	.75 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box
007144	11.5mm	1.5 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box
007145	15.0mm	3.5 cc	7-5 F/2.3-1.7mm	200cm	.035"/.89mm	1/box

■ Metal Stents



GORE® VIABIL® Short Wire Biliary Endoprosthesis

Sterile, Disposable

Recommended Guidewire: .035" (XWire™- 035450SS)

MRI Compatible

Catalog Number	Diameter	Stent Length	Catheter Diameter	O.D.	Unit
Holes					
VSVVH0806	8mm	6cm	200cm	8.5 F	1/box
VSVVH0808	8mm	8cm	200cm	8.5 F	1/box
VSVVH0810	8mm	10cm	200cm	8.5 F	1/box
VSVVH1006	10mm	6cm	200cm	8.5 F	1/box
VSVVH1008	10mm	8cm	200cm	8.5 F	1/box
VSVVH1010	10mm	10cm	200cm	8.5 F	1/box
No Holes					
VSVVN0804	8mm	4cm	200cm	8.5 F	1/box
VSVVN0806	8mm	6cm	200cm	8.5 F	1/box
VSVVN0808	8mm	8cm	200cm	8.5 F	1/box
VSVVN0810	8mm	10cm	200cm	8.5 F	1/box
VSVVN1004	10mm	4cm	200cm	8.5 F	1/box
VSVVN1006	10mm	6cm	200cm	8.5 F	1/box
VSVVN1008	10mm	8cm	200cm	8.5 F	1/box
VSVVN1010	10mm	10cm	200cm	8.5 F	1/box

Warning: The safety and effectiveness of this device for use in the vascular system have not been established.


Gore and Viabil are registered trademarks of W.L. Gore & Associates



SHORT WIRE
BILIARY ENDOPROTHESIS

Intended for palliation of malignant strictures in the biliary tree

■ Metal Stents



GORE® VIABIL® Long Wire Biliary Endoprosthesis

Sterile, Disposable

Recommended Guidewire: .035" (XWire™- 035450SS)


MRI Compatible

Catalog Number	Diameter	Stent Length	Catheter Diameter	O.D.	Unit
Holes					
VH0806200	8mm	60mm	200cm	8.5 F	1/box
VH0808200	8mm	80mm	200cm	8.5 F	1/box
VH0810200	8mm	100mm	200cm	8.5 F	1/box
VH1006200	10mm	60mm	200cm	8.5 F	1/box
VH1008200	10mm	80mm	200cm	8.5 F	1/box
VH1010200	10mm	100mm	200cm	8.5 F	1/box
No Holes					
VN0804200	8mm	40mm	200cm	8.5 F	1/box
VN0806200	8mm	60mm	200cm	8.5 F	1/box
VN0808200	8mm	80mm	200cm	8.5 F	1/box
VN0810200	8mm	100mm	200cm	8.5 F	1/box
VN1004200	10mm	40mm	200cm	8.5 F	1/box
VN1006200	10mm	60mm	200cm	8.5 F	1/box
VN1008200	10mm	80mm	200cm	8.5 F	1/box
VN1010200	10mm	100mm	200cm	8.5 F	1/box



Gore and Viabil are registered trademarks of W.L. Gore & Associates

Metal Stents




Flexxus® Endoscopic Biliary Stent

Sterile, Disposable. Recommended Guidewire: .035" (FXWire™– FX35450S)

MRI Compatible

Catalog Number	Stent		Catheter		Unit
	Diameter	Length	Diameter	Length	
LXB08040	8mm	40mm	7.5 F	190cm	1/box
LXB08060	8mm	60mm	7.5 F	190cm	1/box
LXB08080	8mm	80mm	7.5 F	190cm	1/box
LXB08100	8mm	100mm	7.5 F	190cm	1/box
LXB10040	10mm	40mm	7.5 F	190cm	1/box
LXB10060	10mm	60mm	7.5 F	190cm	1/box
LXB10080	10mm	80mm	7.5 F	190cm	1/box
LXB10100	10mm	100mm	7.5 F	190cm	1/box

■ Fine Needle Aspiration



ClearView™ Endoscopic Ultrasound Needle

Sterile, Disposable

Catalog Number	Description	Gauge	Stylet	Needle Diameter	Units
CV2119R	ClearView™ EUS Needle	19	Round	1.1mm x 2.1mm	5/case
CV1822EB	ClearView™ EUS Needle	22	Extended Bevel	0.7mm x 1.8mm	5/case
CV2722EB	ClearView™ EUS Needle	22	Extended Bevel w/Sheath Stabilizer	0.7mm x 2.7mm	5/case
CV1822R	ClearView™ EUS Needle	22	Round	0.7mm x 1.8mm	5/case
CV2722R	ClearView™ EUS Needle	22	Round w/Sheath Stabilizer	0.7mm x 2.7mm	5/case
CV1825R	ClearView™ EUS Needle	25	Round	0.5mm x 1.8mm	5/case
CV2725R	ClearView™ EUS Needle	25	Round w/Sheath Stabilizer	0.5mm x 2.7mm	5/case

All ClearView™ Needles include lockable aspiration syringe and valve

PET Balloon Dilators




Eliminator® PET Balloon Dilators

Disposable

Catalog Number	Description	Inflated OD		Balloon Length	Catheter Length	Unit
		mm	French			
000340	Esophageal Balloon Dilator	6mm	18 F	8cm	180cm	1/box
000341	Esophageal Balloon Dilator	8mm	24 F	8cm	180cm	1/box
000342	Esophageal Balloon Dilator	10mm	30 F	8cm	180cm	1/box
000343	Esophageal Balloon Dilator	12mm	36 F	8cm	180cm	1/box
000344	Esophageal Balloon Dilator	15mm	45 F	8cm	180cm	1/box
000345	Esophageal Balloon Dilator	18mm	54 F	8cm	180cm	1/box
000380	Esophageal Balloon Dilator	20mm	60 F	8cm	180cm	1/box
000735	Esophageal Advantage Pack (includes 5 ea. of 15mm)			8cm	180cm	5/box
000736	Esophageal Advantage Pack (includes 5 ea. of 18mm)			8cm	180cm	5/box
000379	Esophageal Advantage Pack (includes 5 ea. of 20mm)			8cm	180cm	5/box
000844	Pyloric/Colonic Balloon	6mm	18 F	4cm	230cm	1/box
000845	Pyloric/Colonic Balloon	8mm	24 F	4cm	230cm	1/box
000846	Pyloric/Colonic Balloon	10mm	30 F	4cm	230cm	1/box
000847	Pyloric/Colonic Balloon	12mm	36 F	4cm	230cm	1/box
000848	Pyloric/Colonic Balloon	15mm	45 F	4cm	230cm	1/box
000849	Pyloric/Colonic Balloon	18mm	54 F	4cm	230cm	1/box
000861	Pyloric/Colonic Balloon	20mm	60 F	4cm	230cm	1/box
000858	Pyloric/Colonic Advantage Pack (includes 5 ea. of 15mm)			4cm	230cm	5/box
000859	Pyloric/Colonic Advantage Pack (includes 5 ea. of 18mm)			4cm	230cm	5/box
000867	Pyloric/Colonic Advantage Pack (includes 5 ea. of 20mm)			4cm	230cm	5/box

PET Balloon Dilators



CONMED Balloon Inflation Device

Catalog Number	Description	Unit
100851	CONMED Balloon Inflation Device	1/box

Multi-Stage Inflation Pressures

For Esophageal and Pyloric/Colonic Balloon Dilators

Balloon Size	Stage One	Stage Two	Stage Three*
6mm/18 F	60 psi/4.1 atm	90 psi/6.1 atm	120 psi/8.2 atm
8mm/24 F	60 psi/4.1 atm	90 psi/6.1 atm	120 psi/8.2 atm
10mm/30 F	45 psi/3.1 atm	68 psi/4.0 atm	90 psi/6.1 atm
12mm/36 F	45 psi/3.1 atm	68 psi/4.0 atm	90 psi/6.1 atm
15mm/45 F	30 psi/2.0 atm	45 psi/3.1 atm	60 psi/4.1 atm
18mm/54 F	25 psi/1.7 atm	38 psi/2.6 atm	50 psi/3.4 atm
20mm/60 F	25 psi/1.7 atm	38 psi/2.6 atm	50 psi/3.4 atm

*Maximum recommended inflation pressure

■ Reusable Polyvinyl Dilators

American Dilation System



American Dilation System

Non-Sterile, Reusable

Catalog Number	Description	Size	Unit
000250*	American Dilation System		1/box
000270	Tapered Over-The-Wire Dilator	15 F	1/box
000288	Tapered Over-The-Wire Dilator	18 F	1/box
000271	Tapered Over-The-Wire Dilator	21 F	1/box
000289	Tapered Over-The-Wire Dilator	24 F	1/box
000272	Tapered Over-The-Wire Dilator	27 F	1/box
000273	Tapered Over-The-Wire Dilator	30 F	1/box
000274	Tapered Over-The-Wire Dilator	33 F	1/box
000275	Tapered Over-The-Wire Dilator	36 F	1/box
000290	Tapered Over-The-Wire Dilator	39 F	1/box
000277	Tapered Over-The-Wire Dilator	42 F	1/box
000278	Tapered Over-The-Wire Dilator	45 F	1/box
000279	Tapered Over-The-Wire Dilator	48 F	1/box
000280	Tapered Over-The-Wire Dilator	51 F	1/box
000281	Tapered Over-The-Wire Dilator	54 F	1/box
000291	Tapered Over-The-Wire Dilator	60 F	1/box

*** Each American Dilation System Includes:**

- 15 ea. Tapered over-the-wire dilators
- 1 ea. Reusable Marked Spring Tip Guidewire
- 2 ea. Cleaning adaptor
- 1 ea. Cleaning brush
- 1 ea. Carrying case

■ **CleanGuide™ Disposable Spring Tip Guidewire**

CleanGuide Disposable Spring Tip Guidewire



The NEW Sterile CleanGuide Disposable Marked Spring Tip Guidewire helps reduce the risk of potential cross contamination that could occur from using a reusable guidewire.




Accessories for American Dilation System

DIS150	CleanGuide Disposable Marked Spring Tip Guidewire (Single use)	10/box
000150	Reusable Marked Spring Tip Guidewire	2/box
000249	American Dilator Carrying Case	1/box
000336	American Dilator Cleaning Brush (Single use)	3/box
000337	Cleaning Adaptor	2/box



Carrying Case

Hemostasis Clips




DuraClip™ Repositionable Hemostasis Clip


Sterile, Disposable

Catalog Number	Description	Sheath Length	Min Channel Diameter	Jaw Opening	Qty/Case
DC0165	DuraClip™ for Gastroscopy	165cm	2.8mm	11mm	10 per box
DC0235	DuraClip™ for Colonoscopy	235cm	2.8mm	11mm	10 per box
DC0235W	DuraClip™ for Colonoscopy	235cm	2.8mm	16mm	10 per box

Hemostasis Probes



Three pairs of Positive and Negative Parallel Silver Electrodes




BiCap® Superconductor™ Multi-electrode Bipolar Probes

Sterile, Disposable

Catalog Number	Description	Shaft Size	Shaft Length	Tip Length	Color	Unit
BP-5200A	Bipolar, 6 Silver Electrode, Hemostasis Probe	5 Fr	200cm	6.5mm	Blue	1/box
BP-7300A	Bipolar, 6 Silver Electrode, Hemostasis Probe	7 Fr	300cm	6.5mm	Yellow	1/box
BP-7350A	Bipolar, 6 Silver Electrode, Hemostasis Probe	7 Fr	350cm	6.5mm	Yellow	1/box
BP-7350S	Bipolar, 6 Silver Electrode, Hemostasis Probe	7 Fr	350cm	3.3mm	Yellow	1/box
BP-10300A	Bipolar, 6 Silver Electrode, Hemostasis Probe	10 Fr	300cm	6.5mm	Orange	1/box

■ Hemostasis Probes



Beamer® Argon Snare Probe

Sterile, Disposable

Catalog Number	Description	Unit
Beamsnare-1	Beamer® Snare Probe 2.5mm x 160 (15mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-2	Beamer® Snare Probe 2.5mm x 230 (15mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-3	Beamer® Snare Probe 2.5mm x 230 (30mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-4	Beamer® Snare Probe 2.5mm x 230 (15mm, Oval, Monofilament Firm Loop)	5/case
Beamsnare-5	Beamer® Snare Probe 2.5mm x 230 (30mm, Oval, Monofilament Firm Loop)	5/case



Beamer Argon Probe

Beamer® Argon Probe

Sterile, Disposable

Catalog Number	Description	Unit
A-Beam-1	Beamer® Argon Probe 1.8mm x 160cm	10/case
A-Beam-2	Beamer® Argon Probe 1.8mm x 320cm	10/case
A-Beam-3	Beamer® Argon Probe 2.3mm x 230cm	10/case
A-Beam-4	Beamer® Argon Probe 2.3mm x 320cm	10/case
A-Beam-5	Beamer® Argon Probe 3.2mm x 230cm	10/case



Beamer Argon SideFire Probe

Beamer® Argon SideFire Probe

Sterile, Disposable

Catalog Number	Description	Unit
S-Beam-3	Beamer® Argon SideFire Probe 2.3mm x 230cm	10/case


■ Hemostasis Probes



ABCFlex™ Argon Probes

Sterile, Disposable

Catalog Number	Description	Unit
133023	ABCFlex™ Probe 2.3mm (7Fr.) x 220cm, Original	10/case
133270	ABCFlex™ Probe 2.3mm (7Fr.) x 270cm, Original	10/case
133327	ABCFlex™ Probe 2.3mm (7Fr.) x 270cm, System 7550®	10/case
133000	Extension Adapter for use with ABCFlex™ Probes	1 each



Coaxial Bipolar Adaptors

Non-Sterile, Reusable

Catalog Number	Description	Unit
133001	Coaxial Bipolar Adapter (For use with the System 7500® & System 7550®)	1/box
133002	Coaxial Bipolar Adapter * (For use with the System 5000™ & System 2450™)	1/box

*For use with Everest™ and ValleyLab™ Generators



Active Cord

Non-Sterile, Reusable

Catalog Number	Description	Unit
SN-UAC	Universal Monopolar Active Cord (Fits Olympus style electrode connector)	1/box

■ Disposable Injection Needles



Click-Tip™ Injection Needle

Sterile, Disposable

Catalog Number	Needle Diameter	Needle Length	Catheter OD	Catheter Length	Unit
02-23-230	19 Gauge Injection Needle	4mm	2.3mm	230cm	10
04-23-230	19 Gauge Injection Needle	6mm	2.3mm	230cm	10
03-19-230	22 Gauge Injection Needle	6mm	1.9mm	230cm	10
05-19-230	25 Gauge Injection Needle	4mm	1.9mm	230cm	10
06-19-230	25 Gauge Injection Needle	6mm	1.9mm	230cm	10



Flexitip™ Disposable Sclerotherapy Needles

Sterile, Disposable

Catalog Number	Description	Needle		Sheath		Unit
		Length	Gauge	Length	O.D.	
000215	Optic Yellow Tip Needle	4mm	25 g	160cm	2.3mm	5/box
000216	Optic Yellow Tip Needle	5mm	25 g	160cm	2.3mm	5/box
000218	Optic Yellow Tip Colonoscope Needle	5mm	25 g	230cm	2.3mm	5/box
000198	Standard Needle	4mm	25 g	160cm	2.3mm	5/box
000199	Standard Needle	5mm	25 g	160cm	2.3mm	5/box
000200	Standard Needle	6mm	25 g	160cm	2.3mm	5/box
000138	Standard Colonoscope Needle	5mm	25 g	230cm	2.3mm	5/box

■ Ligation Devices



Stiegmann-Goff™ Endoscopic Esophageal Ligator Kit

Latex

Non-Sterile, Disposable

Catalog Number	Description	Inner Cylinders	Maximum Scope Length	Scope Tip Diameter	Minimum Working Channel O.D.	Unit
200220	Endoscopic Ligator	5	130cm	9 -11mm	2.8mm	5 kits/box
000221	Clearvue™ Endoscopic Ligator	5	130cm	9 -11mm	2.8mm	5 kits/box
000227	Clearvue™ Endoscopic Ligator	10	130cm	9 -11mm	2.8mm	5 kits/box

Each Stiegmann-Goff™ Endoscopic Ligator Includes:

- 1 ea. Monofilament trip wire
- 1 ea. Friction fit adaptors for 9-11mm and 11-13mm scopes
- 5 ea. (200220/000221) inner cylinders with preloaded bands or
- 10 ea. (200221/000227) inner cylinders with preloaded bands



CONMED Endoscopic Overtube

Non-Sterile, Reusable

Catalog Number	Description	O.D.	I.D.	Length	Reuses	Unit
000230*	CONMED Endoscopic Overtube	60 F	45 F	23cm	25	1/box

* For use with the Scope Saver™ Disposable Bite Block, no. 000429

■ Ligation Devices



Stiegmann-Goff™ Bandito™ Hemorrhoid Endoscopic Ligator

Latex

Non-Sterile, Disposable

Catalog Number	Description	Endoscope Compatibility	Maximum Scope Length	Scope Tip Diameter	Minimum Working Channel	Unit
100225	Endoscopic Hemorrhoid Ligator	Colonoscope	210cm	13 -15mm	2.8mm	5 kits/box

- 1 ea. Monofilament trip wire
- 1 ea. Friction fit adaptor for 13-15mm scopes
- 3 ea. Clear inner cylinder with preloaded bands



SpaceBander® Internal Hemorrhoid Band Ligator

Non-Sterile, Disposable

Catalog Number	Description	Unit
HB1020	SpaceBander® Internal Hemorrhoid Band Ligator	20/box

■ Disposable Biopsy Forceps



Oval Cup



Oval Cup with Needle



Alligator Cup



Alligator Cup with Needle

Precisor® EXL™ Coated Disposable Biopsy Forceps

Featuring Channel-Guard™ Low Friction Sheath

Sterile, Coated

Catalog Number	Jaw Configuration	Jaw	Shaft		
		O.D.	Length	Color	Unit
Standard GI					
000381	Alligator Cup	2.3mm	160cm	Yellow	10/box
000382	Alligator Cup with Needle	2.3mm	160cm	Yellow	10/box
000383	Oval Cup	2.3mm	160cm	Yellow	10/box
000384	Oval Cup with Needle	2.3mm	160cm	Yellow	10/box
000385	Alligator Cup	2.3mm	230cm	Blue	10/box
000386	Alligator Cup with Needle	2.3mm	230cm	Blue	10/box
000389	Oval Cup	2.3mm	230cm	Blue	10/box
000388	Oval Cup with Needle	2.3mm	230cm	Blue	10/box
Jumbo GI					
00821C	Oval Cup with Needle	3.1mm	230cm	Blue	10/box
00823C	Alligator Cup with Needle	3.1mm	230cm	Blue	10/box



Oval Cup



Oval Cup with Needle



Alligator Cup



Alligator Cup with Needle

Optibite® Disposable Biopsy Forceps

Sterile, UnCoated

Catalog Number	Jaw Configuration	Jaw Size	Shaft Length	Color	Unit/Box
100571	Alligator cup w/needle	2.3mm	230cm	Blue	10/box
100572	Alligator cup w/o needle	2.3mm	230cm	Blue	10/box
100573	Oval cup w/needle	2.3mm	230cm	Blue	10/box
100574	Oval cup w/o needle	2.3mm	230cm	Blue	10/box
Sterile, Coated					
100575	Alligator cup w/needle	2.3mm	230cm	Blue	10/box
100576	Alligator cup w/o needle	2.3mm	230cm	Blue	10/box
100577	Oval cup w/needle	2.3mm	230cm	Blue	10/box
100578	Oval cup w/o needle	2.3mm	230cm	Blue	10/box

■ Disposable Biopsy Forceps



Oval Cup



Alligator Cup

Pediatric GI

<i>Sterile, Coated</i>		Jaw	Shaft		
Catalog Number	Jaw Configuration	O.D.	Length	Color	Unit
100512	Oval Cup	1.8mm	160cm	Blue	10/box
100513	Alligator Cup	1.8mm	160cm	Blue	10/box



Oval Cup



Alligator Cup

Precisor® Hot Disposable Biopsy Forceps

<i>Sterile, Coated</i>						
Catalog Number	Configuration	Jaw		Shaft		
		O.D.	Length	Coating	Cord	Unit
000852	Oval Cup	2.3mm	230cm	Red	Olympus®	10/box
000854	Alligator Cup	2.3mm	230cm	Red	Olympus	10/box



Oval Cup



Alligator Cup

Precisor® Broncho® Disposable Biopsy Forceps

<i>Sterile, Disposable</i>				
Catalog Number	Description	Shaft Length	Jaw O.D.	Unit
100502	Oval Cup Forceps/Channel-Guard™ Coated Shaft	115cm	1.8mm	10/box
100503	Alligator Cup Forceps/Channel-Guard™ Coated Shaft	115cm	1.8mm	10/box
100506	Oval Cup Forceps	115cm	1.8mm	10/box
100507	Alligator Cup Forceps	115cm	1.8mm	10/box

■ Disposable Percutaneous Biopsy Instrument



Monopty® Percutaneous Biopsy Instrument Kits

Latex

Sterile, Disposable


Catalog Number	Description	Size	Penetration	Notch	Length	Unit
000402*	Monopty® Biopsy Instrument Kit	18 g/1.2mm	22mm	17mm	20cm	5/box
000439*	Monopty® Biopsy Instrument Kt	16 g/1.7mm	22mm	17mm	16cm	5/box
000441*	Monopty® Biopsy Instrument Kit	14 g/2.1mm	22mm	17mm	16cm	5/box

* The syringe plunger in all kits contains latex.

Each Monopty Biopsy Instrument Kit Includes:

- 1 ea. Monopty biopsy instrument includes cm gradations
- 1 ea. 25ml sample vial
- 1 ea. Patient label
- 1 ea. Thumb forceps
- 2 ea. 5.5cm diameter filter paper
- 1 ea. Adhesive bandage
- 1 ea. 24" x 36" fenestrated drape
- 1 ea. Scalpel with No. 11 surgical blade
- 1 ea. 5cc syringe
- 1 ea. 25 gauge x 5/8" needle
- 1 ea. 21 gauge x 1 1/2" needle
- 1 ea. Triple pack povidone-iodine swabs
- 3 ea. 3" x 3" gauze pads
- 1 ea. CSR wrap

■ Disposable GI Cytology Brushes



CONMED Cytology Brushes

Sterile, Disposable

Catalog Number	Description	Scope	Brush Diameter	Sheath O.D.	Length	Unit
000104	Sheathed	Multi-Purpose	1.7mm	1.7mm	160cm	20/box
000108	Sheathed	Pediatric Gastroscope	3.0mm	1.7mm	160cm	20/box
000109	Sheathed	Duodenoscope	3.0mm	1.8mm	200cm	20/box
000110	Sheathed	Colonoscope	3.0mm	2.1mm	230cm	20/box
000196	Sheathed	Adult Gastroscope	3.0mm	2.0mm	160cm	20/box



Cytology Brushes

Sterile, Disposable

Catalog Number	Description	Brush	Sheath		Unit
		Diameter	Length	Diameter	
151	Straight Handle	2.5mm	160cm	1.7mm	20/box
151R	Ring Handle	2.5mm	160cm	1.7mm	20/box
163	Straight Handle	3.0mm	160cm	2.4mm	20/box
158	Straight Handle	3.0mm	230cm	2.4mm	20/box
158R	Ring Handle	3.0mm	230cm	2.4mm	20/box

■ Gastrostomy Feeding Systems



EnTake® PEG System (No Drugs):

Sterile, Disposable

Catalog Number	Description	Technique	French Size	Unit
FESG18ND	Safety	Push	18	1/box
FESG20ND	Safety	Push	20	1/box
FESG24ND	Safety	Push	24	1/box
FESP18ND	Safety	Pull	18	1/box
FESP20ND	Safety	Pull	20	1/box
FESP24ND	Safety	Pull	24	1/box
FEXG20ND	Standard	Push	20	1/box
FEXP20ND	Standard	Pull	20	1/box

Accessories for Gastrostomy Feeding Systems:

Catalog Number	Description	Technique	French Size	Unit
FEDP2400	Dual Port Feeding Adapter	-	24 F	5/box
FEAC1400	Small White Tubing Clamp	-	14 F	10/box
FEAC2000	Large White Tubing Clamp	-	18 F, 20 F and 24 F	10/box

■ Gastrostomy Feeding Systems



EnTake® PEG System

Sterile, Disposable

Catalog Number	Description	Technique	French Size	Unit
FESG1800	Safety	Push	18	1/box
FESG2000	Safety	Push	20	1/box
FESG2400	Safety	Push	24	1/box
FESP1800	Safety	Pull	18	1/box
FESP2000	Safety	Pull	20	1/box
FESP2400	Safety	Pull	24	1/box
FEXG1800	Standard	Push	18	1/box
FEXG2000	Standard	Push	20	1/box
FEXG2400	Standard	Push	24	1/box
FEXP1800	Standard	Pull	18	1/box
FEXP2000	Standard	Pull	20	1/box
FEXP2400	Standard	Pull	24	1/box

EnTake® PEG System includes:

Patient Prep and Access

- 1 ea. Fenestrated drape, 24" x 26"
- 1 ea. Surgical Skin Marker
- 1 ea. Triple pack Chloraprep® Swabsticks**
- 1 ea. Lidocaine, HCl 1% 10cc plastic ampule**
- 1 ea. 5cc syringe*
- 1 ea. Injection Needle
- 1 ea. Scalpel*

- 1 ea. Introducer Needle*
- 1 ea. Snare, disposable grasping

PEG Placement


- 1 ea. Blue cable hoop (Pull Only)
- 1 ea. Guidewire (Push Only)
- 1 ea. PEG Tube (Traction Removal)
- 10 ea. Gauze Pads, 4" x 4"
- 1 ea. Lubricating jelly

- 1 ea. Scissors, blunt end disposable
- 1 ea. Hemostat, disposable curved metal
- 1 ea. Linear bolster
- 1 ea. Star bolster
- 1 ea. Tubing Clamp
- 1 ea. Two Port Feeding adapter
- 1 ea. Gauze pads, 2" x 2" fenestrated
- 1 ea. Ointment, Iodophor PVP**

*Safety features available within Safety PEG System

**Not included in ND kits

■ Disposable Snares




Beamer® Argon Snare Probe

Sterile, Disposable

Catalog Number	Description	Unit
Beamsnare-1	Beamer® Snare Probe 2.5mm x 160 (15mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-2	Beamer® Snare Probe 2.5mm x 230 (15mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-3	Beamer® Snare Probe 2.5mm x 230 (30mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-4	Beamer® Snare Probe 2.5mm x 230 (15mm, Oval, Monofilament Firm Loop)	5/case
Beamsnare-5	Beamer® Snare Probe 2.5mm x 230 (30mm, Oval, Monofilament Firm Loop)	5/case

■ Disposable Snares




Singular™ Polypectomy Snares

Sterile, Disposable

Catalog Number	Loop		Active Cord	Sheath		
	Size	Shape		O.D.	Length	Unit
Firm Loop Wire						
000957	Medium	Crescent	Olympus	2.3mm	230cm	10/box
000959	Medium	Hexagonal	Olympus	2.3mm	230cm	10/box
000963	Large	Oval	Olympus	2.3mm	230cm	10/box
000965	Medium	Oval	Olympus	2.3mm	230cm	10/box
000978	Small (Mini)	Oval	Olympus	2.3mm	230cm	10/box
000991	X-Small	Oval	Olympus	2.3mm	230cm	10/box
Soft Loop Wire						
000977	Large	Oval	Olympus	2.3mm	230cm	10/box
000983	Medium	Oval	Olympus	2.3mm	230cm	10/box
000985	Small (Mini)	Oval	Olympus	2.3mm	230cm	10/box
000987	X-Small	Oval	Olympus	2.3mm	230cm	10/box

■ Disposable Snares




Optimizer® Polypectomy Snares

Sterile, Disposable

Catalog Number	Loop		Active Cord	Sheath		Unit
	Size	Shape		O.D.	Length	
Firm Loop Wire						
000472	Large	Oval	Olympus	2.3mm	230cm	10/box
000474	Medium	Oval	Olympus	2.3mm	230cm	10/box
000476	Small (Mini)	Oval	Olympus	2.3mm	230cm	10/box
000478	X-Small	Oval	Olympus	2.3mm	230cm	10/box
Soft Loop Wire						
000460	Large	Oval	Olympus	2.3mm	230cm	10/box
000462	Medium	Oval	Olympus	2.3mm	230cm	10/box
000464	Small (Mini)	Oval	Olympus	2.3mm	230cm	10/box
000466	X-Small	Oval	Olympus	2.3mm	230cm	10/box

■ Disposable Snares



Orbit-Snare™ Rotatable Polypectomy Snares

Sterile, Disposable

Catalog Number	Loop		Active Cord	Sheath		Wire Type	Unit
	Size	Shape		O.D.	Length		
01-15-230	15mm	Rotatable Oval	Universal	2.5mm	230cm	Firm	10/box
01-25-230	25mm	Rotatable Oval	Universal	2.5mm	230cm	Firm	10/box
01-35-230	35mm	Rotatable Oval	Universal	2.5mm	230cm	Firm	10/box

CONMED Polypectomy Snares

Sterile, Disposable

Catalog Number	Loop		Active Cord	Sheath		Unit
	Size	Shape		O.D.	Length	
Firm Loop Wire						
100600	16mm	Oval	Olympus	2.3mm	230cm	10/box
100601	23mm	Oval	Olympus	2.3mm	230cm	10/box
100602	32mm	Oval	Olympus	2.3mm	230cm	10/box
Soft Loop Wire						
100603	16mm	Oval	Olympus	2.3mm	230cm	10/box
100604	23mm	Oval	Olympus	2.3mm	230cm	10/box
100605	32mm	Oval	Olympus	2.3mm	230cm	10/box




Active Cord

Non-Sterile, Reusable

Catalog Number	Active Cord
SN-UAC	Universal Monopolar Active Cord (Fits Olympus style electrode connector)

■ Specimen Devices



Nakao Spider-Net® Retrieval Devices

Non-Sterile, Disposable

Catalog Number	Description	Net Size	Net Volume	Length	Unit
0230AM	Nakao Spider-Net® Multi-Pack	2.5 x 5.5cm	7cc	230cm	20/box
00230A	Nakao Spider-Net®	2.5 x 5.5cm	7cc	230cm	10/box
0160MG	Nakao Spider-Net® - EGD	2.5 x 5.5cm	7cc	160cm	5/box
0350EN	Nakao Spider-Net® - Enteroscope	2.5 x 5.5cm	7cc	350cm	5/box



Optimizer® Multi-Chamber Polyp Trap

Non-Sterile, Disposable

Catalog Number	Description	Chambers	Unit
CAP132	Multi-Chamber Polyp Trap	4	20/box
CAP132B	Multi-Chamber Polyp Trap	4	120/box

■ Biopsy Valves / Introducer Bite Blocks



ClearGuard™ Biopsy Valves

Sterile, Disposable

Catalog Number	Description	Unit
000522	Biopsy Valve for Olympus and Fujinon	100/case
000523	Biopsy Valve for Pentax	100/case

Biopsy Channel Accessory Introducer

Non-sterile, Disposable

Catalog Number	Description	Unit
000160	Accessory Introducer	6/box



Scope Saver® Bite Blocks

Non-sterile, Disposable

Catalog Number	Description	Lumen Size	Color	Unit
000429 [†]	Large Scope Saver® Disposable Latex Free Bite Block	20 x 27mm	Green	50/box
100429	Small Scope Saver® Disposable Latex Free Bite Block	14 x 20mm	Purple	25/box

[†] Accommodates up to a 60 F American Dilator

CONMED Bite Blocks

Non-sterile, Disposable

Catalog Number	Description	Lumen Size	Color	Unit
001429 [†]	CONMED Large Bite Block Blue Latex Free	20 x 27mm	Blue	100/box

[†] Accommodates up to a 60 F American Dilator

■ Endoscope Valve-Tracking and Pre-Cleaning System


- Helps ensure safe handling of endoscope valves during cleaning, disinfection & storage.
- Provides traceability of valves and other detachable parts, keeping them together with parent endoscope during reprocessing.
- Meets standards of AAMI, American National Standard ST91; 2015 & SGNA Standards of Infection Prevention in Reprocessing of Flexible Gastrointestinal Endoscopes.



CleanCage™ Endoscope Valve-Tracking System

Non-Sterile, Disposable

Catalog Number	Description	Unit
CC8475	Endoscope Valve Cage	100/box



- All-in-one tray, channel cleanser and scope sponge for immediate post-procedure bedside endoscope cleaning.
- Just add water to simplify and standardize your endoscope pre-cleaning processes.
- Meets SGNA Practice Committee Standards of Infection Control.



CleanTray™ Endoscope Pre-Cleaner

Non-Sterile, Disposable

Catalog Number	Description	Unit
CT7357	Bedside Pre-Cleaning Endoscope Kit	100/box

■ CleanGuard™ Disposable Endoscope Valves




CleanGuard™ Disposable Endoscope Valves

Sterile, Disposable

Catalog Number	Description	Unit
200305	Air/Water, Suction and Biopsy Valves for Olympus 140/160/180/190/240/260/290 series Endoscopes	50/case
200310	Air/Water, Suction, Biopsy Valves and Aux Water Jet 140/160/180/190/240/260/290 series Endoscopes	50/case

■ Cleaning Brushes



Channel Master™ Endoscope Cleaning Brushes

Latex **FREE**

Non-Sterile, Disposable - Biopsy Channel Range 2.0-4.2mm, Shaft O.D. 1.8mm

Catalog Number	Scope	Brush		Shaft Length	Unit
		Diameter	Length		
000615	CONMED Double Ended Cleaning Brush	6.4mm	7.4mm	230cm	50/box
000617	CONMED Combination Valve Cleaning Brush	7.9mm/6.4mm	29.2mm/7.4mm	230cm	50/box
000619	CONMED Large Combination Valve Cleaning Brush	9.7mm/6.4mm	45.2mm/7.4mm	230cm	50/box
000639	CONMED Single Ended Cleaning Brush	6.4mm	7.4mm	230cm	50/box



Whistle® Endoscope Cleaning Brush

Non-Sterile, Disposable – Biopsy Channel Minimum 2.8mm, Shaft O.D. 1.8mm

Catalog Number	Scope	Brush		Shaft Length	Unit
		Diameter	Length		
000101	Whistle® Cleaning Brush	6.4mm	13.1mm	230cm	50/box

Contains female luer adaptor, hollow shaft, and port hole for fluid injection to soften channel obstructions.

■ Cleaning Brushes



Blue Bullet® Disposable Cleaning Brushes

Non-Sterile, Disposable

Catalog Number	Scope	Brush		Shaft	
		Diameter	Length	Length	Unit
480	Bronchoscope	3mm	13mm	115cm	50/box
473	Gastro/Duodeno/Colonoscope	3mm	13mm	250cm	50/box
474	Gastro/Duodeno/Colonoscope	4mm	13mm	250cm	50/box
478	Gastro/Duodeno/Colonoscope	5mm	13mm	250cm	50/box
476	Gastro/Duodeno/Colonoscope	6mm	13mm	250cm	50/box
478C	Combination Valve/Channel Cleaning Brush for Gastroscope/Colonoscope	10mm/5mm	38mm/13mm	250cm	50/box
478D	Double End Cleaning Brush for Gastroscope/Colonoscope	5mm	13mm	250cm	50/box

The 400 series Blue Bullet® feature a blue opaque tubing.



Reusable Cleaning Brush

Non-Sterile, Reusable

Catalog Number	Scope	Brush		Shaft	
		Diameter	Length	Length	Unit
380	Bronchoscope	3mm	9.5mm	115cm	1/box
378	Gastro/Duodeno/Colonoscope	6mm	9.5mm	230cm	1/box

The reusable scope cleaning brush features a controllable ring handle and a Teflon® coated spring sheath.

■ Transbronchial Aspiration Needles



WANG® Cytology Needles

Sterile, Disposable

Catalog Number	Biopsy Site	Needle		Unit
		Gauge	Length	
MW-122	Central, Carinal and Mediastinal Regions	22 g	13mm	4/box
MW-222	Central or Peripheral Regions	22 g	13mm	4/box
MW-322	Central, Mediastinal and Hilar Regions	22 g	13mm	4/box
MW-522	Peripheral Region	22 g	13mm	4/box
SW-121	Central, Carinal and Mediastinal Regions	21 g	15mm	4/box
SW-221	Central or Peripheral Regions	21 g	15mm	4/box
SW-521	Peripheral Region	21 g	15mm	4/box
W-220	Central or Peripheral Regions	20 g	15mm	4/box

All WANG® needles have a 1.9mm working diameter and a 140cm working length.



WANG® Histology Needles

Sterile, Disposable

Catalog Number	Biopsy Site	Needle		Unit
		Gauge	Length	
MW-319	Central and Mediastinal Regions	inner - 21 g	3mm	4/box
		outer - 19 g	15mm	
MWF-319	Peripheral Region	19 g	15mm	4/box

All WANG® needles have a 1.9mm working diameter and a 140cm working length.

■ Biopsy Forceps



Oval Cup



Alligator Cup

Precisor Broncho® Disposable Biopsy Forceps

Sterile, Disposable

Catalog Number	Description	Shaft Length	Jaw O.D.	Unit
100502	Oval Cup Forceps/Channel-Guard™ Coated Shaft	115cm	1.8mm	10/box
100503	Alligator Cup Forceps/Channel-Guard™ Coated Shaft	115cm	1.8mm	10/box
100506	Oval Cup Forceps	115cm	1.8mm	10/box
100507	Alligator Cup Forceps	115cm	1.8mm	10/box

Brushes



Cytology Brushes

Sterile, Disposable

Catalog Number	Description	Brush Diameter	Length	Sheath		Unit
				Diameter		
133	Straight Handle	1mm	120cm	1.7mm		20/box
133R	Ring Handle	1mm	120cm	1.7mm		20/box
129	Straight Handle	2mm	120cm	1.7mm		20/box
129R	Ring Handle	2mm	120cm	1.7mm		20/box
149	Straight Handle	3mm	120cm	1.7mm		20/box
149R	Ring Handle	3mm	120cm	1.7mm		20/box



Catalog Number	Description	Scope	Brush Diameter	Length	Sheath		Unit
					Diameter		
000104	Sheathed	Multi-Purpose	1.7mm	160cm	1.7mm		20/box
000107	Sheathed	Bronchoscope	3.0mm	115cm	1.7mm		20/box
000139	Harrell™ (unsheathed)	Bronchoscope	4.0mm	110cm	1.7mm		20/box
000140	Harrell™ (unsheathed)	Bronchoscope	7.0mm	110cm	1.2mm		20/box

■ Brushes



NB-120



PB-120


Needle Tip and Peripheral Cytology Brushes

Sterile, Disposable

Catalog Number	Description	Brush Diameter	Length	Sheath		Unit
				Diameter	Length	
NB-120*	Needle Tip Cytology Brush	1.7mm	11mm	1.9mm	140cm	4/box
PB-120 ^o	Peripheral Radiopaque Tip Cytology Brush	1.7mm	11mm	1.9mm	140cm	4/box

* Needle tip cytology brush is used for retrieving cytology specimens from endobronchial lesions, peripheral lung nodules, or lung masses.

^o The metal tipped sheath of the peripheral brush allows for fluoroscopic visualization during cytologic collection.



Microbiology Brush

Sterile, Disposable

Catalog Number	Description	Brush Diameter	Length	Sheath		Unit
				Diameter	Length	
130	Microbiology Brush	1mm	110cm	1.8mm		5/box

Designed to obtain contamination free specimens, the brush is protected behind a wax plug which is expelled before sampling.

■ Cleaning Brushes



Blue Bullet® Disposable Cleaning Brush

Non-Sterile, Disposable

Catalog Number	Scope	Brush		Length	Shaft Length	Unit
		Diameter				
480	Bronchoscope	3mm		13mm	115cm	50/box


The 400 series Blue Bullet® features a blue opaque tubing.



Channel Master™ Endoscope Cleaning Brush

Non-Sterile, Disposable

Catalog Number	Scope	Biopsy Channel Range	Shaft Length	Unit
000610	Bronchoscope	2.0 - 4.2mm	160cm	50/box *



Reusable Cleaning Brush

Non-Sterile, Reusable

Catalog Number	Scope	Brush		Shaft Length	Unit
		Diameter	Length		
380	Bronchoscope	3mm	9.5mm	115cm	1/box

The reusable scope cleaning brush features a controllable ring handle and a Teflon® coated spring sheath.

Generators



Beamer® System CE600

Catalog Number	Description
CE200	ESU
CB200	Argon Module
CC200	Cart
CE200-A05	Beamer® Double Foot Pedal with Swap-Mode Button and Foot Lift Bar
CB200-A02	Argon Pressure Reducer
CB200-A01	Beamer® Connecting Cable for Argon Probes
CE200-A02	Beamer® Single Foot Pedal with Lift Bar
CE200-A04	Connecting Cable for Reusable Return Electrodes
CC200-A10	Irrigation Pump Mounting Shelf
SN-UAC	Universal Monopolar Active Cord

Generators



BiCap® III ElectroSurgical Generator

Catalog Number	Description
60-8512-SYS	BiCap® III System with Monopolar and Bipolar Footswitches*
60-8512-120	BiCap® III ElectroSurgical Generator
60-6700-001	Monopolar Footswitch
60-5103-001	BiPolar Footswitch
60-8510-001	BiCap® III Pump Mounting Bracket

*Includes 60-8512-120, 60-6700-001 and 60-5103-001.




Irrigation Pump

Catalog Number	Description	Unit
CETEGP-100	Irrigation Pump for Endoscopy Fits on Beamer™ Cart with Beamer® System CE600 or as a Stand Alone Irrigation System	1/box

Each CETEGP-100 includes:

1 ea. Foot Switch, 1 ea. Power Cord, 1 ea. Bottle Holder, 1 ea. Pump Head

■ Argon Probes



Beamer® Argon Snare Probes

Sterile, Disposable

Catalog Number	Description	Unit
Beamsnare-1	Beamer® Snare Probe 2.5mm x 160 (15mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-2	Beamer® Snare Probe 2.5mm x 230 (15mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-3	Beamer® Snare Probe 2.5mm x 230 (30mm, Firm Spiral Braided Oval Loop)	5/case
Beamsnare-4	Beamer® Snare Probe 2.5mm x 230 (15mm, Oval, Monofilament Firm Loop)	5/case
Beamsnare-5	Beamer® Snare Probe 2.5mm x 230 (30mm, Oval, Monofilament Firm Loop)	5/case



Beamer® Argon Probes

Sterile, Disposable

Catalog Number	Description	Unit
A-Beam-1	Beamer® Probe 1.8mm x 160cm	10/case
A-Beam-2	Beamer® Probe 1.8mm x 320cm	10/case
A-Beam-3	Beamer® Probe 2.3mm x 230cm	10/case
A-Beam-4	Beamer® Probe 2.3mm x 320cm	10/case
A-Beam-5	Beamer® Probe 3.2mm x 230cm	10/case



Beamer® Argon SideFire Probe

Sterile, Disposable

Catalog Number	Description	Unit
S-Beam-3	Beamer® Argon SideFire Probe 2.3mm x 230cm	10/case



ABCFlex™ Argon Probes

Sterile, Disposable

Catalog Number	Description	Unit
133023	ABCFlex™ Probe 2.3mm (7Fr.) x 220cm, Original	10/case
133270	ABCFlex™ Probe 2.3mm (7Fr.) x 270cm, Original	10/case
133327	ABCFlex™ Probe 2.3mm (7Fr.) x 270cm, System 7550®	10/case
133000	Extension Adapter for use with ABCFlex™ Probes	1 each

Please be advised that the following devices or components provided by CONMED contain natural rubber latex, which may cause allergic reactions in some individuals:

Product	Catalog Number	Device/Component
Stiegmann-Goff™ Endoscopic Ligator Kit	000221, 000227, 200220-221	"O" ring contains latex
Stiegmann-Goff™ Bandito™ Endoscopic Ligator	100225	"O" ring contains latex
DuraGlide3™ Stone Removal Balloon	007123-007125, 007143-007145	Latex balloon
DuraGlide™ Stone Removal Balloon	000635, 000637-638, 000646-647	Latex balloon
Disposable Percutaneous Biopsy Instruments	000402, 000439, 000441	Syringe Plunger

Catalog Number	Page	Catalog Number	Page	Catalog Number	Page
000101.....	37	000336.....	16	000639.....	37
000104.....	26	000337.....	16	000735.....	13
000104.....	41	000340.....	13	000736.....	13
000107.....	41	000341.....	13	000844.....	13
000108.....	26	000342.....	13	000845.....	13
000109.....	26	000343.....	13	000846.....	13
000110.....	26	000344.....	13	000847.....	13
000138.....	20	000345.....	13	000848.....	13
000139.....	41	000347.....	7	000849.....	13
000140.....	41	000348.....	7	000852.....	24
000150.....	16	000349.....	7	000854.....	24
000160.....	34	000379.....	13	000858.....	13
000196.....	26	000380.....	13	000859.....	13
000198.....	20	000381.....	23	000861.....	13
000199.....	20	000382.....	23	000867.....	13
000200.....	20	000383.....	23	000957.....	30
000215.....	20	000384.....	23	000959.....	30
000216.....	20	000385.....	23	000963.....	30
000218.....	20	000386.....	23	000965.....	30
000221.....	21	000388.....	23	000977.....	30
000227.....	21	000389.....	23	000978.....	30
000230.....	21	000402.....	25	000983.....	30
000249.....	16	000429.....	24	000985.....	30
000250.....	15	000439.....	25	000987.....	30
000270.....	15	000441.....	25	000991.....	30
000271.....	15	000460.....	31	00230A.....	33
000272.....	15	000462.....	31	001429.....	34
000273.....	15	000464.....	31	007104.....	3
000274.....	15	000466.....	31	007105.....	3
000275.....	15	000472.....	31	007123.....	8
000277.....	15	000474.....	31	007124.....	8
000278.....	15	000476.....	31	007125.....	8
000279.....	15	000478.....	31	007143.....	8
000280.....	15	000522.....	34	007144.....	8
000281.....	15	000523.....	34	007145.....	8
000288.....	15	000610.....	43	00821C.....	23
000289.....	15	000615.....	37	00823C.....	23
000290.....	15	000617.....	37	01-15-230.....	32
000291.....	15	000619.....	37	01-25-230.....	32

CATALOG NUMBER INDEX

Catalog Number	Page	Catalog Number	Page	Catalog Number	Page
01-35-230.....	32	100577.....	23	478D.....	38
02-23-230.....	20	100578.....	23	480.....	38
0160MG.....	33	100600.....	32	60-5103-001.....	45
0230AM.....	33	100601.....	32	60-6700-001.....	45
025260S.....	6	100602.....	32	60-8510-001.....	45
025450A.....	6	100603.....	32	60-8512-120.....	45
025450S.....	6	100604.....	32	60-8512-SYS.....	45
03-19-230.....	20	100605.....	32	7104AC.....	3
0350EN.....	33	100851.....	14	7105AC.....	3
035260A.....	6	129.....	41	A25261S.....	5
035260AS.....	6	129R.....	41	A25481S.....	5
035260S.....	6	130.....	42	A25261A.....	5
035260SS.....	6	133.....	41	A25481A.....	5
035450A.....	6	133000.....	19	A-Beam-1.....	18
035450AS.....	6	133000.....	46	A-Beam-1.....	46
035450S.....	6	133001.....	19	A-Beam-2.....	18
035450SS.....	6	133002.....	19	A-Beam-2.....	46
04-23-230.....	20	133023.....	19	A-Beam-3.....	18
05-19-230.....	20	133270.....	19	A-Beam-3.....	46
050906.....	4	133270.....	46	A-Beam-4.....	18
050907.....	4	133327.....	19	A-Beam-4.....	46
06-19-230.....	22	133327.....	46	A-Beam-5.....	18
100225.....	22	133R.....	41	A-Beam-5.....	46
100429.....	34	149.....	41	Beamsnare-1.....	18
100502.....	24	149R.....	41	Beamsnare-1.....	29
100502.....	40	151.....	26	Beamsnare-2.....	18
100503.....	24	151R.....	26	Beamsnare-2.....	29
100503.....	40	158.....	26	Beamsnare-3.....	18
100506.....	24	158R.....	26	Beamsnare-3.....	29
100506.....	40	163.....	26	Beamsnare-4.....	18
100507.....	24	200220.....	21	Beamsnare-4.....	29
100507.....	40	200305.....	36	Beamsnare-5.....	18
100512.....	24	200310.....	36	Beamsnare-5.....	29
100513.....	24	378.....	38	BP-10300A.....	17
100571.....	23	380.....	38	BP-5200A.....	17
100572.....	23	473.....	38	BP-7300A.....	17
100573.....	23	474.....	38	BP-7350A.....	17
100574.....	23	476.....	38	BP-7350S.....	17
100575.....	23	478.....	38	CAP132.....	33
100576.....	23	478C.....	38	CAP132B.....	33

Catalog Number	Page	Catalog Number	Page	Catalog Number	Page
CB200	44	FESP24ND	27	T2043F	2
CB200-A01	44	FEXG1800	28	T3043F	2
CB200-A02	44	FEXG18ND	27	VH0806200	10
CC200	44	FEXG2000	28	VH0808200	10
CC200-A10	44	FEXG20ND	27	VH0810200	10
CC8475	35	FEXG2400	28	VH1006200	10
CE200	44	FEXG24ND	27	VH1008200	10
CE200-A02	44	FEXP1800	28	VH1010200	10
CE200-A04	44	FEXP2000	28	VN0804200	10
CE200-A05	44	FEXP20ND	27	VN0806200	10
CETEGP-100	45	FEXP2400	28	VN0808200	10
CT7357	35	FX35260A	5	VN0810200	10
CV1822EB	12	FX35260S	5	VN1004200	10
CV1822R	12	FX35450A	5	VN1006200	10
CV1825R	12	FX35450S	5	VN1008200	10
CV2119R	12	HB1020	22	VN1010200	10
CV2722EB	12	LXB08040	11	VSWVN0804	9
CV2722R	12	LXB08060	11	VSWVN0806	9
CV2725R	12	LXB08080	11	VSWVN0808	9
DC0165	17	LXB08100	11	VSWVN0810	9
DC0235	17	LXB10040	11	VSWVN1004	9
DC0235W	17	LXB10060	11	VSWVN1006	9
DIS150	16	LXB10080	11	VSWVN1008	9
F922SW	8	LXB10100	11	VSWVN1010	9
F922LW	8	MW-122	39	VSWVH0806	9
FEAC1400	27	MW-222	39	VSWVH0808	9
FEAC2000	27	MW-319	39	VSWVH0810	9
FEDP2400	27	MW-322	39	VSWVH1006	9
FESG1800	28	MW-522	39	VSWVH1008	9
FESG18ND	27	MWF-319	39	VSWVH1010	9
FESG2000	28	NB-120	42	W-220	39
FESG20ND	27	PB-120	42		
FESG2400	28	S-BEAM-3	18		
FESG24ND	27	SN-UAC	19		
FESP1800	28	SN-UAC	32		
FESP18ND	27	SN-UAC	46		
FESP2000	28	SW-121	39		
FESP20ND	27	SW-221	39		
FESP2400	28	SW-521	39		



Gastroenterology and Pulmonary Product Catalog

525 French Road, Utica, NY 13502

Customer Service: 1-866-4CONMED

International Sales: +1-727-214-3000

Fax: 1-800-438-3051

Website: www.CONMED.com

Domestic Email: customerexperience@conmed.com

International Email: internationalorders@conmed.com