



Knee Preservation System™

Multiple continuous stitches for repair of any size meniscal tear.

Sequent® Meniscal Repair Device



- 6 Stitch
- 3 Stitch
- 2 Stitch




MULTIPLE
CONTINUOUS STITCHES



ALL-INSIDE
STAY-INSIDE



KNOTLESS
TECHNIQUE



To learn more about these and other innovative products, call 800-448-6506 or visit ConMed.com.

KNEE
PRESERVATION
SYSTEM™



The Sequent® Meniscal Repair Device is an All-Inside, Stay-Inside system comprised of PEEK implants and Hi-Fi® suture that can create up to six continuous stitches for a knotless repair in any configuration.


This technology provides surgeons a meniscal repair technique with the opportunity for a faster repair. It also reduces meniscal trauma and the risk of chondral damage.

Target any size meniscal tear pathology, including oblique, longitudinal, transverse, horizontal and bucket handle with ease, using as few as two or as many as six continuous stitches.

Sequent® devices come in three different size configurations, allowing surgeons to choose the right tool to repair virtually any meniscal tear using the easy-to-master PRDT™ technique.

PRDT™: *Pierce, Rotate, Deploy, Tension*
TECHNIQUE

Ultimate Fixation Strength



*Data on File

KNEE
PRESERVATION
SYSTEM™

As the chart demonstrates, *Ultimate Fixation Strength* performance of the Sequent® Meniscal Repair Device surpasses similar products in the field. This is a clear clinical benefit for patient outcomes.

Sequent® Meniscal Repair Device

Sequent Meniscal Repair Device, Curved Needle, **7 Implants**

MR007C

Sequent Meniscal Repair Device, Curved Needle, **4 Implants**

MR004C

Sequent Meniscal Repair Device, Curved Needle, **3 Implants**

MR003C

Sequent Disposable Kit (Suture Cutter, Entry Cannula)

SCEC047

Sequent Disposable Suture Cutter

SC047D

Sequent Reusable Entry Cannula

EC047



ConMed Linvatec
Knee Preservation System™
Sequent® Meniscal Repair Device

525 French Road
Utica, New York 13502

Local 727-392-6464
Toll Free 800-237-0169

ConMed.com
customer_service@conmed.com



Sequent® Anchors are comprised of PEEK-OPTIMA® polymer from Invivo® Biomaterial Solutions. PEEK-OPTIMA and INVIVIO are registered trademarks of Invivo Ltd. All Rights are Reserved.